

♣ Color code by flow rate

No.	Flow rate	Color code	Color
1	0.5ml/hr	Pantone 1788C	0.5ml/hr
2	1.0ml/hr	Pantone 3005C	1.0ml/hr
3	1.5ml/hr	Pantone 298C	1.5ml/hr
4	2.0ml/hr	Pantone Yellow C	2.0ml/hr
5	2.5ml/hr	Pantone 4695C	2.5ml/hr
6	3.0ml/hr	Pantone 5473C	3.0ml/hr

No.	Flow rate	Color code	Color
7	4.0ml/hr	Pantone 165C	4.0ml/hr
8	5.0ml/hr	Pantone 354C	5.0ml/hr
9	6.0ml/hr	Pantone 572C	6.0ml/hr
10	8.0ml/hr	Pantone Cool gray 4C	8.0ml/hr
11	10.0ml/hr	Pantone 267C	10.0ml/hr